

AIDAPT BATHROOMS LTD.

Lancots Lane
Sutton Oak
St. Helens
WA9 3EX



Solo Bed Transfer Aid (With Strap)

The Solo Bed Transfer Aid assists transfer into and out of one side of the bed, without affecting other occupants. The hand grip sits above the bed and can be adjusted in height to suit, while the base frame is large enough to provide stability.

Features:

- Height adjustable grip with safety rail
- Large base for improved stability
- High quality easy to clean epoxy resin finish
- 5 height adjustment settings
- Tubular steel frame

PRODUCT SPECIFICATION

Depth (mm).....	680
Height (mm).....	580-680
Width (mm)	600
Product Dimensions (mm)	580x600x680
Net weight (kg)	3.31
Colour	White
Configuration.....	Transfer Aid & Strap set
Handle Width (mm).....	325
Maximum User Weight (kg)	127
Maximum User Weight (stones).....	20
Overall Depth (mm).....	680
Width of Base (mm).....	600



Product Code: VY438

Barcode: 5021196598220

Registered Office:
AIDAPT BATHROOMS LTD.
LANCOTS LANE, SUTTON OAK, ST. HELENS, WA9 3EX